1 Command Line Dictionary

 $yes|no\{no\}$ acf automask low-pass limit(in Å) amsklp angular error(in degrees) angerr append $yes|no\{no\}$ astigmatism error(in degrees) astigerr yes|no avg bfactor bfac bfactor error bfacerr binarize{no} bin image size(in pixels) box boxfiles.txt boxtab $yes|no\{no\}$ center chunksz chunk size cluster2extract class $clip2box{256}$ clip yes|homo|no|{no} clustvalid $shc_clustering_nclsX.txt$ comlindoc compare $yes|no{no}$ corner corner size(in pixels) $yes|no\{no\}$ countvox cs spherical aberration constant(in mm) $\{2.7\}$ yes|no|flip|mul{no} ctf ctfsq $yes|no\{no\}$ $yes|no\{no\}$ ctfstats cube side(in pixels) cube deferr defocus error(in microns) defocus defocus(in microns) text file with CTF info deftab density dens directory dir discrete $yes|no\{no\}$ $yes|no{yes}$ diversify list of oritabs doclist dynlp $yes|no{yes}$ 1st Euler angle e1 e2 2nd Euler angle 3d Euler angle e3 edge size(in pixels) edge $yes|no{yes}$ eo $yes|no\{no\}$ errify $yes|no\{no\}$ even body of file fbody movies.txt filetab find Fourier index filename fname fraction frac fraction of amplitude contrast $\{0.07\}$ fraca from particle index fromp

 $\begin{array}{ll} \texttt{fsc} & & \texttt{fsc_state1.bin} \\ \texttt{ft2img} & & \texttt{yes}|\texttt{no}\{\texttt{no}\} \\ \texttt{guinier} & & \texttt{yes}|\texttt{no}\{\texttt{no}\} \end{array}$

hfun sigm|tanh|lin{sigm}

hist var2plot

hp high-pass limit(in Å)

inner inner mask radius(in pixels)
kv acceleration voltage{300}
label class|state|subclass{class}
lp low-pass limit(in Å){20}
lpstart starting low-pass limit
lpstop stay at this low-pass limit

masscen $yes|no\{no\}$

maxits maximum number of iterations {500}

minp minimum cluster population

mirr $no[2d]3d\{no\}$

moldiam molecular diameter(in Å)
msk mask radius(in pixels)

 $\begin{array}{ll} {\tt mskfile} & {\tt msk.ext} \\ {\tt msktype} & {\tt hard|soft\{soft\}} \end{array}$

mul shift multiplication factor mw molecular weight(in kDa)

ncls nr of clusters
ndiscrete nr of discrete oris
ndocs nr of documents

neg yes|no

newbox new box size(in pixles)
nframes nr of movie frames

 $\begin{array}{lll} \mbox{noise} & yes|no\{no\} \\ \mbox{noise_norm} & yes|no \\ \mbox{norm} & yes|no\{no\} \end{array}$

nptcls nr of particles

nran size of random sample

nspace nr of orientations in discrete space

nstates nr of conformational states

nthr nr of threads (CPUs within a socket)

 $\begin{array}{lll} \text{nvox} & \text{nr of voxels} \\ \text{odd} & \text{yes}|\text{no}\{\text{no}\} \\ \text{oritab} & \text{orientations} \\ \text{oritab2} & \text{orientations} \\ \text{outfile} & \text{output text file} \\ \text{outstk} & \text{output image stack} \\ \end{array}$

 $\begin{array}{lll} \text{outvol.ext} & \text{outvol.ext} \\ \text{pgrp} & \text{cn}|\text{dn}|t|\text{o}|\text{i}\{\text{c1}\} \\ \text{phrand} & \text{yes}|\text{no}\{\text{no}\} \\ \text{plot} & \text{yes}|\text{no}\{\text{no}\} \end{array}$

pspecsz box size of power spectrum

refine no|shc{no}

refs initial references.ext

 $\begin{array}{lll} \text{rnd} & \text{yes}|\text{no}\{\text{no}\}\\ \text{round} & \text{yes}|\text{no}\{\text{no}\}\\ \text{scale} & \text{scale factor}\\ \text{shalgn} & \text{yes}|\text{no}\{\text{no}\}\\ \end{array}$

 $\begin{array}{ll} \texttt{sherr} & \text{shift error(in pixels)} \\ \texttt{smpd} & \text{sampling distance(in \mathring{A})} \end{array}$

snr signal2noise ratio

soften $yes|no\{no\}$

split nr of partitions to split into

srch_inpl yes|no{yes}
startit start iteration nr
state state index
stats yes|no|print{no}
stk input particle stack
stk2 input particle stack nr 2
stk3 input particle stack nr 3

thres threshold

top stop particle index

trs origin shift range[-trs,trs](in pixels)

 $\begin{array}{ll} \text{vis} & \text{yes}|\text{no}\{\text{no}\}\\ \text{vol1} & \text{invol.ext}\\ \text{vol2} & \text{invol2.ext} \end{array}$

which_iter iteration number width pixels width